

CRImage, a package for classifying cells and caculating tumour cellularity

Henrik Failmezger

April 21, 2010

1 Load the package

The package is loaded by the following command:

```
> library(CRImage)
```

See the functions of the package EBIImage to read, write and manipulate an image.

2 Threshold an image

To create a binary image of a coloured image the function `calculateThreshold` can be used. This function calculates Otsu threshold on an image.

```
> f = system.file("data", "exImg2.jpg", package = "CRImage")
> exImgB = readImage(f)
> indexWhitePixel = which(exImgB[, , 1] > 0.85 & exImgB[, , 2] >
+ 0.85 & exImgB[, , 3] > 0.85)
> exImgB = channel(exImgB, "gray")
> t = calculateThreshold(as.vector(exImgB[-indexWhitePixel]))
> exImgB[which(exImgB >= t)] = 1
> exImgB[which(exImgB < t)] = 0
```

The image is read with `readImage`. `indexWhitePixel` finds white pixels in the image. These pixels are considered as background and excluded from thresholding. `channel(exImgB, "gray")` the image is converted to grayscale and thresholded with `calculateThreshold`. The binary image is created by finding pixels with gray values larger or smaller than the threshold and colouring them black or white.

3 Segment an image

An image can be segmented to find cells in the image, using the function `segmentImage`.




Figure 1: Thresholded image.

```
> f = system.file("data", "exImg2.jpg", package = "CRImage")
> segmentationValues = segmentImage(f)
```

The image is converted to grayscale and thresholded. Morphological opening is used to delete clutter and to smooth the shapes of the cells. The watershed algorithm is used to separate clustered cells. `segmentationValues[[1]]` holds the original image, `segmentationValues[[2]]` holds the segmented image, `segmentationValues[[3]]` holds features, which were calculated for the segmented objects. The segmented objects can be drawn by the function `display(paintObjects(segmentedIm` (see `EBImage`).

4 Creating a training set

To classify cells in images first an appropriate training set has to be created.

```
f = system.file("data", "exImg.jpg", package="CRImage")
trainingValues=createTrainingSet(filename=f)
```

The list `trainingValues` returns two values. The first value is an image which can be displayed by:

```
display(trainingValues[[1]])
```

is an image in which every segmented cell is numbered. The second value is an table with features for every cell:

| index | class | g.x | g.y | g.s | g.p | g.pdm | g.pdsd | g.effr... |
|-------|-------|----------|-----------|-----|-----|-----------|-----------|--------------|
| 1 | <NA> | 148.2203 | 4.855932 | 118 | 36 | 5.721049 | 1.0130550 | 6.128668... |
| 2 | <NA> | 160.2719 | 4.763158 | 114 | 38 | 5.659780 | 1.0331399 | 6.023896... |
| 3 | <NA> | 183.7975 | 3.101266 | 79 | 35 | 4.593585 | 0.9575905 | 5.014627... |
| 4 | <NA> | 196.3500 | 4.242857 | 140 | 43 | 6.424591 | 1.4433233 | 6.675581... |
| 5 | <NA> | 271.5694 | 2.680556 | 72 | 29 | 4.504704 | 1.2490794 | 4.787307... |
| 6 | <NA> | 338.5221 | 6.530973 | 113 | 35 | 5.601531 | 1.0849651 | 5.997418... |
| 7 | <NA> | 456.0179 | 2.946429 | 112 | 39 | 5.726556 | 1.8101368 | 5.970821... |
| 8 | <NA> | 556.3778 | 9.018519 | 270 | 81 | 9.196894 | 2.4687870 | 9.270581... |
| 9 | <NA> | 575.1777 | 10.935950 | 484 | 101 | 11.959995 | 2.1378649 | 12.412171... |
| 10 | <NA> | 592.7391 | 1.724638 | 69 | 41 | 5.145094 | 2.2782539 | 4.686511... |

To create the training set class values for the cells have to be inserted in the column `class`:

| index | class | g.x | g.y | g.s | g.p | g.pdm | g.pdsd | g.effr |
|-------|-----------|----------|----------|-----|-----|----------|-----------|----------|
| 1 | normal | 148.2203 | 4.855932 | 118 | 36 | 5.721049 | 1.0130550 | 6.128668 |
| 2 | malignant | 160.2719 | 4.763158 | 114 | 38 | 5.659780 | 1.0331399 | 6.023896 |

The values in the column `index` match the numbers for the cells in the image. You can save the table as tab separated by:

```
write.table(trainingData,file="path",sep="\t",rownames=F)
```




Figure 2: SegmentedImage.

and open it for example in a spreadsheet program, however be careful not to shift the columns. You do not have to assign a class value to every cell but do ensure that there are enough examples for every class. Class values can be numbers (1,2,3) or strings (normal, malignant). You can choose at most 10 classes. If you want to use the kernel smoothing approach for classification you can only choose two classes.

5 Creating the classifier

The command:

```
> f = system.file("data", "trainingData.txt", package = "CRImage")
> trainingData = read.table(f, header = TRUE)
> classifierValues = createClassifier(trainingData, topo = FALSE)
> classifier = classifierValues[[1]]
```

creates the classifier. The classifier is a Support Vector Machine provided by the package e1017.

6 Classification of cells

After having created the classifier, cells in another image can be classified.

```
> f = system.file("data", "exImg2.jpg", package = "CRImage")
> classValues = classifyCells(classifier, filename = f, KS = TRUE)
```

7 Calculation of cellularity

If tumour images are processed, the cellularity of the tumour can be calculated. The image is first segmented and the cell types are classified. Afterwards, the cellularity of the tumour is determined. The function needs to know, which class value should be the tumour class.

```
> t = system.file("data", "trainingData.txt", package = "CRImage")
> trainingData = read.table(t, header = TRUE)
> classifier = createClassifier(trainingData, topo = FALSE)[[1]]
> f = system.file("data", "exImg2.jpg", package = "CRImage")
> exImg = readImage(f)
> cellularityValues = calculateCellularity(f, classifier = classifier,
+   cancerIdentifier = "1")
```




Figure 3: Heatmap of cellularity values.




Figure 4: Classified Image.